




**LL65600 Shoulder Straps** - lift straps incl. carbiners and shoulder cushions - intended to relieve and facilitate manual lifting. The hooks are of the best quality, used by mountaineers. Protection is available to wrap around the hooks for a softer grip.



#### How to use the Shoulder Straps

Adjust the length of the straps by moving the hooks up or down in the loops on the straps. You can have a straight back or bend backwards depending on which loops you choose. Test to find the best position. It should feel good and secure.

The hooks are also positioned according to the height of the people doing the lifting. A short and a tall person put the hooks at different heights. If lifting without the protection, VikOm, the seatlift handles next to the hook are held.

**VikOm gives a nicer grip.**

For longer lifts, such as down to the water or going up stairs, the Shoulder Straps can be a sensible complement to the liftseat. They fit just as well on the original Lätta Lyft as on the bathing model.

More information is available to download under the respective headings at [www.lattalyft.se](http://www.lattalyft.se) or Facebook Lätta Lyft

Lätta Lyft / Svart på Vitt AB  
+46 707 525820

Guarantee: 2 years from the date of the invoice.

Check that the shoulder straps do not have any defects.

If damage does occur, the straps must be replaced. This product is CE marked.

Manufactured in Sweden

**Svart på Vitt Reklam och Design AB / Lätta Lyft made in Sweden**

Box 4067 S- 443 12 Gråbo Sweden

Tel +46 707 525820 [www.lattalyft.se](http://www.lattalyft.se) [www.facebook.com/lattalyft](https://www.facebook.com/lattalyft)